



Gluteus Maximus Injection Simulator

IM training for the three ways of injection

This simulator is designed to practice for the three ways of injection into the hip and the buttocks. The incorporated indicator shows the accuracy of needle placement and its depth

Product Supervision

Dr. Hiroshi Koyama,
Medical Director, General Internal Medicine
National Hospital Organization Kyoto Medical
Center



FEATURES


- 1 | The life-sized lower torso manikin consist of three layers: skin, muscular and skeleton
- 2 | Greater trochanter, anterior superior iliac spine and scrum can be palpated to identify the injection sites
- 3 | Three methods for intramuscular injection: Quarter-Trichotomy, Hochstetter Method and Clark Method
- 4 | The indicator shows the accuracy of the needle placement and its depth
- 5 | Correct needle placement can be confirmed by being able to inject liquid
- 6 | No needle trace remains

SKILLS

| Intramuscular injection



Hochstetter Method



Clark / Sachtleben Method

DESCRIPTIONS

SPECIFICATIONS

Size: W32 x D35 x H20 cm Weight: 3.2kg / 7.05lbs
W12.6 x D14 x H8 in.

SET INCLUDES

- 1 lower torso manikin
- 1 stand for Hochstetter method
- 1 indicator
- 1 set of sensor pads
- 1 instruction manual

MATERIALS

Soft resin
Latex free

REPLACEMENT PARTS

- 11213-010 skin
- 11213-020 sensor for puncture pad
- 11213-030 sensor for Hochstetter's site
- 11213-040 sensor for Clark's site

RECOMMENDED DEVICES

IM needle 21-23G