

M86

11324-010

Pediatric Patient Care Model (5-6 years old)



Life-size model of 5-6-year-old female child

Ideal training partner for pediatric nursing care

FEATURES

- 1 | Natural articulation of joints
- 2 | A variety of possible patient positioning
- 3 | Maintenance of the patient position
- 4 | Seamless, one-piece waterproof structure
- 5 | Soft and life-like touch of the skin

SKILLS

- | Immobilization technique
- | Fixation of a butterfly needle
- | Changing of clothes of a patients with connected catheters
- | Oxygen tent
- | Body position changing
- | Changing clothes, Bed bathing

DESCRIPTIONS

SPECIFICATIONS

Size: 98 cm / 38.6 in.
Head circumference: 33 cm / 13 in.

Weight: 7.7 kg / 16.9 lbs

SET INCLUDES

- 1 | pediatric manikin
- 1 | pajama
- 1 | instruction manual

MATERIALS

Soft resin, metal internal frame
Latex free

M73

11324-000

Pediatric Patient Care Model (7-10 months)



Life-size model of 7-10-month-old female baby

FEATURES

- 1 | Natural articulation of joints
- 2 | A variety of possible patient positioning
- 3 | Seamless, one-piece waterproof structure
- 4 | Soft and life-like touch of the skin

SKILLS

- | Immobilization technique
- | Fixation of a butterfly needle
- | Changing of clothes of a patients with connected catheters
- | Oxygen tent
- | Body position changing
- | Changing clothes, Bed bathing

DESCRIPTIONS

SPECIFICATIONS

Size: 70 cm / 27.5 in.
Head circumference: 33 cm / 13 in.

Weight: 5.6 kg / 12.5 lbs

SET INCLUDES

- 1 | pediatric manikin
- 1 | pajama
- 1 | instruction manual

MATERIALS

Soft resin, metal internal frame
Latex free