

Pediatric Patient Care Simulator

A comprehensive day-to-day, long-term care manikin for pediatric patients

Provides a variety of features for practicing essential skills and simulating pediatric patient care scenarios for nurses and caregivers

Product Supervision

Medical Association of Osaka Prefecture
Pediatric Medical Care Exploratory Committee








FEATURES

- 1 | Life-like 6 years old model
- 2 | Able to sit in an upright position without support
- 3 | Oral, nasal and tracheal airway suction interchangeable female and male genitalia for catheterization training

SKILLS

- | Positioning
- | Patient handling
- | Airway suction:
 - oral, nasal or through a tracheostomy opening
- | Percutaneous Endoscopic Gastrostomy (PEG) care
- | Feeding tube management
- | Male & female catheterization

IMAGES



DESCRIPTIONS

SPECIFICATIONS

Size: 125cm / 49.2 in. Weight: 13kg / 28.6lbs

SET INCLUDES

1 pediatric manikin	1 strap-on IV venipuncture trainer
1 head & neck anatomy plate	1 uniform (gym clothes)
1 male & female genitalia	1 set of catheters (sample)
1 syringe	1 lubricant
1 irrigation bag	1 instruction manual

MATERIALS

Soft resin, hard resin
Latex free

REPLACEMENT PARTS

11355-020	head and mask explanation panel
11355-030	valve for catheterization
11355-040	replacement mouth/ nose part

* Instruments are not included.